








THUASNE®

Rebel® Lock

fr Genouillère ligamentaire rigide articulée avec mécanisme de blocage.....	4
en Hinged rigid ligament knee brace with locking mechanism.....	4
de Starre Knieführung orthese mit Arretiermechanismus.....	5
nl Rigide scharnierende ligamentaire kniebrace met vergrendelingsmechanisme.....	6
it Ginocchiera per legamenti rigida articolata con meccanismo di blocco.....	7
es Rodillera rígida articulada para ligamentos con mecanismo de bloqueo.....	7
pt Joelheira ligamentar rígida articulada com mecanismo de bloqueio.....	8
da Stift artikuleret knæbind til ledbånd med låsemekanisme.....	9
fi Nivelletty polven jäykkä nivelsidetuki, jossa on lukitusmekanismi.....	9
sv Ledat styvt knäskydd för ledband med blockeringsmekanism.....	10
el Επιγονατίδα συνδέσμων άκαμπτη αρθρωτή με μηχανισμό ασφάλισης.....	11
cs Pevná kolenní vazová ortéza s kloubem se zajišťovacím mechanismem.....	12
pl Sztywna, przegubowa orteza stawu kolanowego z mechanizmem blokującym.....	12
lv Stingra ceļa locītavas saišu ortoze ar eņģēm un bloķēšanas mehānismu.....	13
lt Standus raiščių antkelis su jungtimis ir blokavimo mechanizmu.....	14
et Liigenduvate küljetugevduste ja lukustusmehhanismiga põlvetugi.....	14
sl Toga opornica za kolenske vezi s sklepi in zapornim mehanizmom.....	15
sk Pevná kolenná ortéza na podporu väzov so systémom uzamykania.....	16
hu Ízületes merev ligament térdrögzítő reteszelő mechanizmussal.....	16
bg Твърда съчленена накованка със заключващ механизъм.....	17
ro Genunchieră ligamentară rigidă articulată cu mecanism de blocare.....	18
ru Жесткий лигаментарный ортез для коленного сустава с шарнирами и механизмом блокировки.....	19
hr Čvrsta zglobna ortoza za ligamente koljena s mehanizmom za blokiranje.....	19
zh 带有锁定结构的硬质关节韧带护膝.....	20
ar داعم ركبة للأربطة صلب ومفصلي مزود بآلية للحصر.....	21

						
(a)	(b)	(c)				
cm	in	cm	in	cm	in	
XS	31 - 39	12 - 15	7,5 - 9	3 - 3.5	28 - 32	11 - 12.5
SM	39 - 47	15 - 18.5	9 - 10	3.5 - 4	31 - 35	12 - 13.5
MD	47 - 53	18.5 - 21	10 - 11,5	4 - 4.5	34 - 38	13 - 15
LG	53 - 60	21 - 23.5	11,5 - 12,5	4.5 - 5	37 - 40	14.5 - 15.5
XL	60 - 64	23.5 - 25	12,5 - 14	5 - 5.5	38 - 43	15 - 17

fr	Stabilisation	Immobilisation
en	Stabilisation	Immobilisation
de	Stabilisierung	Ruhigstellung
nl	Stabilisatie	Immobilisatie
it	Stabilizzazione	Immobilizzazione
es	Estabilización	Inmovilización
pt	Estabilização	Imobilização
da	Stabilisering	Støtte
fi	Stabilioi	Immobilisoi
sv	Stabilisering	Immobilisering
el	Σταθεροποίηση	Ακινητοποίηση
cs	Stabilizace	Znehybnění
pl	Stabilizacja	Unieruchomienie
lv	Stabilizēšana	Imobilizācija
lt	Stabilizavimas	Imobilizavimas
et	Stabiliseerimine	Liikumatuks muutmine
sl	Stabiliziranje	Imobilizacija
sk	Stabilizácia	Imobilizácia
hu	Stabilizálás	Rögzítés
bg	Стабилизация	Обездвижване
ro	Stabilizare	Imobilizare
ru	Стабилизация	Иммобилизация
hr	Stabilizacija	Imobilizacija
zh	稳定	固定
ar	الاستقرار	التثبيت

fr	(a) Circonférence de la cuisse (15 cm au-dessus du genou)	(b) Largeur du genou	(c) Circonférence du mollet (15 cm au-dessous du genou)
en	(a) Thigh circumference (15 cm above the knee)	(b) Knee width	(c) Calf circumference (15 cm below the knee)
de	(a) Oberschenkelumfang (15 cm über dem Knie)	(b) Kniebreite	(c) Wadenumfang (15 cm unter dem Knie)
nl	(a) Omvang van de dij (15 cm boven de knie)	(b) Breedte van de knie	(c) Omvang van de kuit (15 cm onder de knie)
it	(a) Circonferenza della coscia (15 cm al di sopra del ginocchio)	(b) Larghezza del ginocchio	(c) Circonferenza del polpaccio (15 cm al di sotto del ginocchio)
es	(a) Circunferencia del muslo (15 cm por encima de la rodilla)	(b) Ancho de la rodilla	(c) Circunferencia de la pantorrilla (15 cm por debajo de la rodilla)
pt	(a) Circunferência da coxa (15 cm acima do joelho)	(b) Largura do joelho	(c) Circunferência da barriga da perna (15 cm abaixo do joelho)
da	(a) Lårets omkreds (15 cm over knæet)	(b) Knæets bredde	(c) Læggens omkreds (15 cm under knæet)
fi	(a) Reiden ympäry (15 cm polven yläpuolelta)	(b) Polven leveys	(c) Pohkeen ympäry (15 cm polven alapuolelta)
sv	(a) Omkrets runt låret (15 cm ovanför knäet)	(b) Knäets bredd	(c) Omkrets runt vad (15 cm under knäet)
el	(a) Περίμετρος του μηρού (15 εκ. πάνω από το γόνατο)	(b) Περίμετρος του γόνατος	(c) Περίμετρος της κνήμης (15 cm κάτω από το γόνατο)
cs	(a) Obvod stehna (15 cm nad kolennem)	(b) Obvod kolena	(c) Obvod lýtky (15 cm pod kolennem)
pl	(a) Obwód uda (15 cm powyżej kolana)	(b) Szerokość kolana	(c) Obwód łydki (15 cm poniżej kolana)
lv	(a) Augšstilba apkārtmērs (15 cm virs ceļgala)	(b) Ceļgala platums	(c) Apakšstilba apkārtmērs (15 cm zem ceļgala)
lt	(a) Šlaunies apimtis (15 cm virš kelio)	(b) Kelio plotis	(c) Blauzdos apimtis (15 cm žemiau kelio)
et	(a) Reie ümbermõõt (15 cm põlvest kõrgemal)	(b) Põlve ümbermõõt	(c) Säãremarja ümbermõõt (15 cm põlvest allpool)
sl	(a) Obseg stegna (15 cm nad kolenom)	(b) Širina kolena	(c) Obseg goleni (15 cm pod kolenom)
sk	(a) Obvod stehna (15 cm nad kolenom)	(b) Šírka kolena	(c) Obvod lýtky (15 cm pod kolenom)
hu	(a) Comb körfogata (a térd felett 15 cm-rel)	(b) Térd szélessége	(c) Vádli körfogata (a térd alatt 15 cm-rel)
bg	(a) Обиколка на бедрото (15 cm над коляното)	(b) Ширина на коляното	(c) Обиколка на прасеца (15 cm под коляното)
ro	(a) Circumferința coapsei (la 15 cm deasupra genunchiului)	(b) Lățimea genunchiului	(c) Circumferința gambei (la 15 cm sub genunchi)
ru	(a) Диаметр окружности бедра (на 15 см выше колена)	(b) Ширина колена	(c) Обхват голени (на 15 см ниже колена)
hr	(a) Opseg bedra (15 cm iznad koljena)	(b) Širina koljena	(c) Opseg lista (15 cm ispod koljena)
zh	(a) 大腿周长 (膝盖上方15厘米处)	(b) 膝关节宽度	(c) 小腿周长 (膝盖下方15厘米处)
ar	(a) محيط الفخذ (15 سم فوق الركبة)	(b) عرض الركبة	(c) محيط ريلة الساق (15 سم تحت الركبة)

MD

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Conservser cette notice.

en HINGED RIGID LIGAMENT KNEE BRACE WITH LOCKING MECHANISM

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Product available in left or right version.

Available in 5 sizes.

Rebel' Lock is available in length 35,5 cm/14 in.

The device is composed of :

- a frame ensuring the brace's rigidity,
- the TM5+ hinge ➊ mimicking the knee's natural movement, protected by a cover,
- 2 condylar pads ➋,
- 4 posterior half-straps (Ⓐ Ⓢ Ⓣ Ⓦ) (including a synergistic suspension strap Ⓐ), and 2 anterior half-straps (ⓈⓉⓌ).

A kit with 4 additional condylar pads ➋ (2 thicknesses: 1cm/¼in and 1,3 cm/½ in), 2 Allen keys (Ⓢ and Ⓣ), 2 flexion washers ➌ and 2 flexion screws ➍ is included in the brace's box.

Depending on country of sale, additional accessories/spare parts could be available.

Composition

Textile components: elastane - polyamide - ethylene vinyl acetate - silicone. Rigid components: aluminium - stainless steel - polyoxymethylene - vinyl - polypropylene - polyamide - epoxy resin - brass - polyurethane.

Properties/Mode of action

Stabilisation of knee joint ligaments provided by the rigid hinged side uprights.

The brace can be locked at 0° and 5° of extension.

Flexion can be adjusted at 30°, 60° and 90°.

The brace is maintained on the leg thanks to:

- the Synergistic Suspension Strap ②,
- the shape of the tibial bolster (rounded lateral part, flat medial part) that prevents rotation,
- the Compression/Suspension pack providing improved hold on the thigh.

Indications

Severe joint instability/laxity (including for knee osteoarthritis).
Conservative treatment of knee ligament injuries and/or ruptures (cruciate and/or lateral ligaments).
Post-operative immobilisation or/and rehabilitation.
Post-traumatic immobilisation.

Contraindications

Do not use the product if the diagnosis has not been confirmed.
Do not use in the event of known allergy to any of the components.
Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.
History of venous or lymphatic disorders.
Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.
Do not use the product if it is damaged.
Choose the appropriate size to fit the patient, referring to the size chart.
It is recommended that a healthcare professional supervises the first application.
Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.
For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional. In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging machine.
Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.
Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

Contouring the brace:

To make the knee brace more comfortable and improve its hold, the uprights can be shaped.

Apply a slight pressure in the desired direction.
Repeat the process on the other upright if needed.
If subsequent shaping of the uprights is needed, consult a healthcare professional.

The use of thicker condylar pads, provided with the brace, may enhance the fit and provide more proprioceptive pressure to the sides of the knee.

Adjusting the strap length:

If one or more of the knee brace straps is too long, each strap can be cut to the desired length.

To do this, remove the self-fastening hook tab from the end of the strap, cut the strap to the desired length, and reattach the self-fastening hook tab to the end of the strap ①.

Take care to not cut any strap too short.
Additionally, any comfort pad that is attached to the inside of a strap must be removed and cut shorter if the pad interferes with your ability to completely tighten the strap.

Adjusting the flexion limitation ①:

Flexion can be limited to 30°, 60° and 90°.
To change this setting, please follow the following instructions, **to be repeated with the same limitation on both hinges**.

To set flexion at 30°, 60° or 90°, remove the the condylar pad to access the flexion stop settings on the inside of the hinge.

Place a flexion screw ④ through the flexion washer ③ and insert both into the desired flexion setting.

Use the smallest Allen key ⑤ to securely tighten the screw.

Adjusting the extension locking ② ⑥:

Extension can be locked to 0° and 5°.
The brace's default locking position is 0°.
To change this setting, please follow the following instructions, **to be repeated with the same limitation on both hinges**.

To change locking position from 0° to 5°, use the largest Allen key ⑤ to remove both screws 3A.

Slide thigh shell backwards until both holes align at the "5" setting.
Reinsert and tighten the screws.

To unlock the hinge pull up on the catch 3B until it clicks ②.
Hinge will move freely between extension and flexion settings.
To lock the hinge, press button (b) at the front of lock 3B and straighten leg completely.

Hinge will lock in the chosen extension position (0° or 5°).
The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.

Both hinges MUST be adjusted to the same angle. There can be damage to the hinges, and compromises (including injury) to the patient, if the settings are not the same on both hinges.

Positioning the brace on the leg:

- Have the patient sit on the edge of a chair, and have the patient bend his/her knee to a position of approximately 30-60° of flexion.
Ensure that the four straps that cross the back of the leg are open (② ③ ④ ⑤), and to ease application, fold the self-fastening straps back.

- Place the knee brace on the leg.
Ensure that the brace size is appropriate by checking that there is no space between the brace and the leg. If there is any gapping, or if the brace is too tight, if possible, conform the frame of the brace or retake measurements of the leg and select the appropriate size (see size table).

- Identifying knee center:
The condylar pads attached to the inside of the hinges should press against the sides of your knee.

The center of the hinges against the sides of the knee must be aligned with the upper third of the patella.

NOTE: It is always better to put your brace on a little too high than too low.
- First tighten and fasten the synergistic suspension strap ② located above the calf. It is the most important strap for preventing the brace from slipping down your leg.

By tightening the synergistic suspension strap ② above the calf muscle in the flexion fold behind the knee you will ensure that the hinges are at the correct level on the leg.

- Next fasten the anterior tibial strap ③ but not too tight.
The anterior-posterior position of the hinges can be modulated by adjusting straps ③ and ④.

The hinges and uprights should be aligned with the midline of the sides of the leg, or just posterior to midline (2nd third of the leg in the sagittal plane).
After you pre-set the anterior strap, instruct the patient to leave it closed when putting on or taking off the brace.

- Next fasten the lower strap ⑤, and then the two straps behind the thigh ⑥ and ⑦, and finish by fastening the anterior thigh strap ⑧.
Strap chafes can be angled to ensure the straps fit flat against the back side of the leg.

- As a final step, push down on the brace to make sure it is strapped securely and stays in place.

We also recommend to stretch out the leg and retighten the lower thigh strap ⑥ to achieve maximum control.
If the patient feels the need to tighten the straps during an activity, make sure to tighten them following the steps indicated in the fitting instructions.

Care/Maintenance

Close the self-fastening tabs before washing. Hand wash. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Hinges ①:

The hinges on the brace are pre-lubricated in the factory.
If sand, dirt or water gets inside the hinges, they may require lubricating again. If you notice the hinges not gliding smoothly, a few drops of a synthetic lubricant can be applied.

Wipe off any excess lubricant before wearing the brace to prevent stains on clothing.

Straps:
After prolonged use, if the fibers on your strap do not adhere as well to the self-fastening tab, cut the strap shorter so the self-fastening tab adheres to a section of the strap that has fresher fibers.

If this is not possible, you should contact the medical provider who fit your brace.

Pads:

The brace is lined with padding that provides a comfortable interface between the leg and the shells.
The straps also have pads.
Do not remove these pads from the brace or straps.

Wipe the pads after each use to remove any moisture and let the pads air dry. You can also clean the pads with a mild anti-bacterial soap and rinse them off with fresh water.

Do not wash the pads in a machine and do not dry in a tumble dryer.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS
Thuasne offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:

- six months for the straps, padding and condylar pads,
- one year for the rigid components.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:

- misuse of the product or any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,
- damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne entity.

Any warranty claim will first be reviewed by Thuasne to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

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©Thuasne - 2047001 (2023-06)