
**THUASNE®**

# LYMPHATREX Revolution

<b>fr</b>	Dispositif de compression sur mesure patient.....	3
<b>en</b>	Made-to-measure patient compression device .....	5
<b>de</b>	Für Patienten maßgefertigtes Kompressionsprodukt .....	7
<b>nl</b>	Compressiehulpmiddel op maat van de patiënt .....	9
<b>it</b>	Dispositivo di compressione su misura paziente .....	12
<b>es</b>	Dispositivo de compresión a medida del paciente .....	14
<b>pt</b>	Dispositivo de compressão adaptado ao paciente .....	16
<b>da</b>	Kompressionsudstyr efter patientens mål.....	19
<b>fi</b>	Potilaan mittojen mukaan valmistettu kompressiotuote.....	21
<b>sv</b>	Individanpassad kompressionsprodukt .....	23
<b>el</b>	Τεχνολογικό προϊόν συμπίεσης προσαρμοσμένο στον ασθενή.....	25
<b>cs</b>	Kompresivní pomůcka na míru .....	28
<b>pl</b>	Wyrób uciskowy na miarę.....	30
<b>lv</b>	Pacientam pielāgots kompresijas izstrādājums.....	32
<b>lt</b>	Pagal paciento užsakymą pagaminta kompresinė priemonė .....	34
<b>et</b>	Patsiendile eritellimusel valmistatav kohandatud kompressioontoode.....	37
<b>sl</b>	Kompresijski pripomoček po meri pacienta.....	39
<b>sk</b>	Kompresívna pomôcka vyrobená na mieru pre pacienta .....	41
<b>hu</b>	A páciens méretére készült kompressziós eszköz.....	43
<b>bg</b>	Компресивно изделие, адаптирано за пациента.....	45
<b>ro</b>	Dispozitiv de compresie realizat la comandă pentru pacient.....	48
<b>ru</b>	Компрессионные изделия по меркам пациента .....	50
<b>hr</b>	Kompresijski proizvod po mjeri bolesnika.....	53
<b>zh</b>	患者定制按压装置.....	55
<b>ja</b>	カスタムメイドの圧迫装置.....	56
<b>ko</b>	환자 맞춤형 압박스타킹.....	58
<b>ar</b>	تجهيزه ضغط على مقياس المريض.....	63

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instruções de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - 使用方法 - 사용설명서 - إرشادات الاستعمال

fr

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Conserver cette notice.

en

## MADE-TO-MEASURE PATIENT COMPRESSION DEVICE

### Description/Destination

The Lymphatrex Revolution range is a range of flat knit compression devices made to the measurements of the patient's limb.

Use: The device is intended for use on the lower or upper limbs of adult or pediatric patients suffering from at least one of the indications listed below, directly on the skin.

This range is available in several models:

Upper limb:	Lower limb:
<ul style="list-style-type: none"> <li>• Armsleeve</li> </ul>	<ul style="list-style-type: none"> <li>• Below-knee</li> <li>• Thigh-high</li> <li>• Leg sleeve</li> </ul>

For measurement, refer to the measurement sheets available on the Thuasne website.

Different finishing options are available according to your healthcare professional's recommendations.

### Composition

Polyamide, elastane, cotton, polyester.

Anti-slip system: silicone, elastane, polyamide, cotton, polyester.

Compression pad option: ethylene-vinyl acetate, polyamide, elastane.

### Properties/Mode of action

Lymphatrex Revolution devices provide compression, applying pressure to the limb to promote lymphatic circulation, help prevent the formation of edema and maintain its volume. They are available in multiple compression classes according to different referentials:

ASQUAL (France)	RAL (Germany)
<ul style="list-style-type: none"> <li>• class 2 (15 - 20 mmHg)</li> <li>• class 3 (20 - 36 mmHg)</li> <li>• class 4 (&gt; 36 mmHg)</li> </ul>	<ul style="list-style-type: none"> <li>• class I (18 - 21 mmHg)</li> <li>• class II (23 - 32 mmHg)</li> <li>• class III (34 - 46 mmHg)</li> </ul>

Refer to the product label to know the compression standard and the associated compression class.

### Indications

- Treatment of acute or chronic edema (e.g: primary or secondary lymphedema of the upper or lower limb).
- Venous diseases of the lower limbs (the recommended compression classes according to French Health Authority (HAS) are indicated in brackets):
  - Heavy legs, restless legs, evening edema (15 - 20 mmHg).
  - Telangiectasia or reticular veins (< 3 mm) (15 - 20 mmHg).
  - Following sclerotherapy for telangiectasia or reticular veins (< 3 mm) (15 - 20 mmHg).
  - Prevention of venous thrombosis in pregnant women or post-partum (15 - 20 mmHg).
  - Varicose veins (≥ 3 mm) (15 - 20 mmHg and 20 - 36 mmHg).
  - Following sclerotherapy or surgery for varicose veins (≥ 3 mm) (15 - 20 mmHg and 20 - 36 mmHg).

- Chronic venous edema (20 - 36 mmHg).
  - Trophic disorders: pigmentation or eczema, chronic lipodermatosclerosis or chronic white atrophy (20 - 36 mmHg).
  - Healed venous ulcer (20 - 36 mmHg and > 36 mmHg).
  - Active venous ulcer (> 36 mmHg).
  - Chronic venous disease associated with pregnancy or the post-partum period (20 - 36 mmHg).
  - Prevention and/or treatment of venous thrombosis (lower limb) and post-thrombotic syndrome (20-36 mmHg).
- Symptomatic treatment of lipedema.

### Contraindications

Do not use in case of major dermatological conditions on an affected limb.

Do not use in the event of known allergy to any of the components.

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of septic thrombosis.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of the affected limb.

#### Specific to the lower limbs:

- Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6.
- Do not use in the event of decompensated heart failure.
- Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg).
- Do not use in case of extra-anatomical bypass of the affected limb.

#### Specific to the upper limbs:

- Do not use in case of brachial plexus condition.
- Do not use in case of vasculitis of the extremities.

### Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

In the event of a venous ulcer, a Doppler examination is necessary before using the device in order to ensure there is no associated severe arterial disease.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional.

Do not apply the product directly to an open wound without a dressing. Take care not to move the dressing or create folds in it when applying or removing the device.

These compression devices are designed for daytime use.

Due to the constant evolution of the pathology, regular follow-up with a healthcare professional is recommended to check that the compression device is still suitable for the pathology.

For some pathologies (or situations) such as venous thrombosis, the device must be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

In the event of lipedema, compression can be used to help relieve the symptoms. Seek advice from a healthcare professional for compression and other solutions for managing lipedema.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

- Do not use when swimming.
- Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product.

Take care not to damage the device when fitting it, particularly with your nails.

Check that the anti-slip system can be applied to healthy, unbroken skin, if applicable.

- For hygiene, security and performance reasons, do not re-use the product for another patient.

### Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

### Instructions for use/Application

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Lymphatrex Revolution compression devices should be worn directly on the skin.

Ask a health professional for advice.

The adequate pressure obtention may need the superposition/combination with other compression devices.

#### Putting on the armsleeve:

- Insert the hand into the garment **1** and turn the garment back on itself up to the opening provided for the hand. **2**
- Insert the hand into the garment then slide the garment carefully along the limb up to the armpit. **3**
- Avoid pulling on the garment more than is necessary.

- Adjust the garment such that no folds remain along the limb.
- Once the garment is fully in place, massage the limb to ensure even compression. ④

Putting on the below-knee, leg sleeves and thigh-high:

- Insert the hand into the garment and turn it back on itself up to the heel. ⑤ (not applicable for the leg sleeve)
- Insert the foot into the garment until the heel is correctly positioned. ⑥ (not applicable for the leg sleeve)
- Slide the garment carefully along the limb to below the knee (below-knee) or to below the gluteal fold (thigh-high). ⑦
- Avoid pulling on the garment more than is necessary.
- Adjust the garment such that no folds remain along the limb.
- Below-knee or thigh-high with anti-slip system: when putting the garment on, do not pull on the anti-slip part to avoid damaging it.
- Once the garment is fully in place, massage the limb to ensure even compression. ⑧

**Care/Maintenance**

- We recommend washing your garment before the first use.
- We recommend washing your garment every day, or after each use.
- Lymphatrex Revolution compression devices are designed for daily use over a 6-month period (i.e. 90 washings).
- After this time, garments should be replaced to ensure the treatment is effective and technical performance is not affected.

Machine washing:

- Machine washable at 30°C (delicate program).
- Use a laundry detergent without fabric softener. Do not use aggressive products (e.g.: chlorinated products).
- Put the garment in a washing net to prevent other garments (especially those with hook-and-loop fasteners or zippers) from damaging the knitted fabric of your orthosis.
- Wring out without twisting.
- Dry away from any direct heat source (radiator, sun, etc.).
- Do not tumble-dry.
- Do not iron.

**Storage**

Store at room temperature, out of direct sunlight.

**Disposal**

Dispose of in accordance with local regulations.

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©Thuasne - 2053901 (2024-05)
**THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186