THUASNE®

Ortel C4 Rigid

fr	Collier cervical rigide.....	3
en	Rigid cervical collar.....	3
de	Feste Halskrause.....	3
nl	Stijve nekbrace.....	3
it	Collare cervicale rigido.....	3
es	Collarín cervical rígido.....	4
pt	Colar cervical rígido.....	4
da	Stiv nakkekrave.....	4
fi	Jäykkä kaulatuki.....	4
sv	Rigid nackkrage.....	5
el	Αυχενικό κολάρο σκληρό.....	5
cs	Pevný krční límec.....	5
pl	Kołnierz sztywny do unieruchomienia szyi.....	5
lv	Stingra kakla ortoze.....	6
lt	Standus kaklo įtvaras.....	6
et	Jäik kaelatugi.....	6
sl	Opornica za vrat, toga.....	6
sk	Pevný krčný golier.....	7
hu	Kemény nyakmerevítő.....	7
bg	Твърда цервикална яка.....	7
ro	Guler cervical rigid.....	7
ru	Жесткий шейный воротник.....	8
hr	Tvrda ortoza za vratnu kralježnicu.....	8
zh	说明/用途.....	8
ar	طوق للعنق صلب.....	8

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Odnoście użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

S	7,50 - 8,75 cm
M	8,75 - 10 cm
L	10 - 11,25 cm
XL	11,25 - 12,50 cm

© Thuasne - 6166104 (2023-03)
© Studio Caterin

THUASNE
120, rue Marius Auphan
92300 Levallois-Perret
France

www.thuasne.com

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RIGID CERVICAL COLLAR

Description/Destination

2 anatomically-shaped moulded shells, made of hypoallergenic foam:
- 1 anterior shell, with chin and sternal support. ①
- 1 posterior shell, with occipital and upper dorsal support. ②
Anterior tracheal opening. ③
Precise self-fastening adjustment. ④
X-ray compatible. ⑤
The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Polyamide - polyethylene - polypropylene - polyester - acrylic.

Properties/Mode of action

Immobilisation of the cervical spine.

Indications

Cervical radiculopathy.
Cervical spondylosis.
Stable cervical fractures.
Postoperative use.
Whiplash.

Contraindications

Do not apply the product in direct contact with broken skin.
Do not use in the event of known allergy to any of the components.
Do not use in the event of traumatic injuries: severe sprain, unstable fracture.
Do not use for prolonged wear without medical supervision.
Do not use the product if the diagnosis has not been confirmed.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

In the event of discomfort, significant hindrance, pain, abnormal sensations, consult a healthcare professional.
Choose the appropriate size to fit the patient, referring to the size chart.
Verify the product's integrity before every use.
It is recommended that a healthcare professional supervises the first application.
This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.
For hygiene and performance reasons, do not re-use the product for another patient.
Do not use the device if it is damaged.
Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).
It is recommended to adequately tighten the device to achieve support without excessive compression.
Rigid cervical collars should be removed as soon as clearance of the cervical spine is achieved, even in intubated patients or obtunded patients (with reduced level of consciousness).
Do not wear the device when driving a vehicle.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters,etc.) or wounds of various degrees of severity.
Wearing a cervical collar may cause some weakness of the neck muscles.
Wearing a cervical collar may raise intracranial pressure.
Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Apply the back part of the collar over the patient's neck. ④
Apply the front part of the collar with the chin resting on the rest provided for this purpose.
Center the collar to ensure alignment in a neutral position. ⑥ ⑥
Close the self-fastener: the front part overlaps the back part. ⑥
Ensure the strap is not fastened too tight. ⑥

Care/Maintenance

See product label. Close the self-fastening tabs before washing.
Hand wash. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine).
Do not tumble-dry. Do not iron. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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