







**THUASNE®**

## Biflex Self Adjust

en	Adjustable compression device for arm and hand.....	6
nl	Verstelbaar compressiehulpmiddel voor arm en hand.....	6
sv	Justerbar kompressionsanordning för arm och hand.....	7
pl	Regulowana opaska uciskowa na ramię lub dłoń.....	8
fr	Dispositif compressif ajustable pour bras et main.....	9
de	Verstellbare Kompressionsbandage für Arm und Hand.....	10
it	Dispositivo di compressione regolabile per braccio e mano.....	11
es	Dispositivo compresivo ajustable para brazo y mano.....	12
pt	Dispositivo de compressão ajustável para o braço e a mão.....	13
da	Justerbar kompressionsanordning til arm og hånd.....	14
fi	Säädettävä kompressiotuote käsivarteen ja käteen.....	15
no	Justerbar kompresjonsinnretning for arm og hånd.....	15
el	Ρυθμιζόμενο σύστημα συμπίεσης βραχίονα και άκρου χεριού.....	16
cs	Nastavitelná kompresní pomůcka pro paži a ruku.....	17
lv	Pielāgojams kompresijas izstrādājums rokai un plaukstai.....	18
lt	Reguliuojama kompresinė priemonė rankai ir plaštakai.....	19
et	Reguleeritav kompressiooniseade käsivarrele ja käele.....	20
sl	Nastavljivi kompresijski pripomoček za roko in dlan.....	21
sk	Nastavitelná kompresívna pomôcka na rameno a ruku.....	21
hu	Állítható kompressziós eszköz karra és kézre.....	22
bg	Регулируемо компресиращо изделие за ръка (рамо до китка и китка до пръсти).....	23
ro	Dispozitiv de compresie ajustabil pentru braț sau mână.....	24
hr	Podesiva oprema za kompresiju ruke i šake.....	25
zh	可调节的手臂或手部压力装置.....	26
ja	腕と手用の調節可能な圧迫装置.....	27
ko	과 손을 위한 조절 가능한 압박 장치.....	27
ar	تجهيزه ضاغطة قابلة للتعديل للذراع واليد.....	29







Arm - Arm - Arm - Ramię - Bras - Arm - Braccio - Brazo - Braço - Arm - Käsivarsi - Arm -  
 - Βραχίονας - Raže - Roka - Ranka - Käsivars - Roka - Rameno - Кар - Рѣка (рамо до  
 китка) - Braç - Ruka - 手臂 - 腕 - 팔 - للذراع

	1	15 - 19 cm 6 - 7.5"	22 - 28 cm 9 - 11"	25 - 32 cm 10 - 12.5"
	2	16 - 20 cm 6.5 - 8"	24 - 30 cm 9.5 - 12"	29 - 36 cm 11.5 - 14"
	3	17 - 21 cm 6.5 - 8"	26 - 32 cm 10 - 12.5"	33 - 40 cm 13 - 15.5"
	4	17 - 21 cm 6.5 - 8"	28 - 34 cm 11 - 13.5"	37 - 44 cm 14.5 - 17.5"
	5	18 - 22 cm 7 - 8.5"	30 - 36 cm 12 - 14"	41 - 48 cm 16 - 19"
	6	19 - 23 cm 7.5 - 9"	31 - 37 cm 12.5 - 14.5"	45 - 52 cm 18 - 20.5"

	cG	cD	cC
<b>en</b>	Axillary arm circumference	Forearm circumference (widest point)	Wrist circumference
<b>nl</b>	Omtrek van de arm onder de oksel	Omtrek van de onderarm (grootste omtrek)	Omtrek van de pols
<b>sv</b>	Armens omfång vid axillen (G)	Omkrets underarm (på det bredaste stället)	Omkrets handled
<b>pl</b>	Obwód ramienia pod pachą	Obwód przedramienia (w najszerszym miejscu)	Obwód nadgarstka
<b>fr</b>	Circonférence de bras axillaire	Circonférence de l'avant-bras (au plus fort)	Circonférence de poignet
<b>de</b>	Armumfang auf Achselhöhe	Unterarmumfang (an der stärksten Stelle)	Handgelenkumfang
<b>it</b>	Circonferenza braccio a livello ascellare	Circonferenza avambraccio (nel punto più ampio)	Circonferenza polso
<b>es</b>	Contorno axilar	Contorno del antebrazo (parte más gruesa)	Contorno de la muñeca
<b>pt</b>	Circunferência axilar do braço	Circunferência do antebraço (na parte maior)	Circunferência do pulso
<b>da</b>	Overarmens omkreds ved armhulen	Underarmens omkreds (det kraftigste sted)	Håndleddets omkreds
<b>fi</b>	Olkavarren ympärysmitta kainalon kohdalta	Kyynärvarren ympärysmitta (paksuimmasta kohdasta)	Ranteen ympärysmitta
<b>no</b>	Omkrets ved armhulen	Forarmens omkrets (på det tykkeste)	Håndleddets omkrets
<b>el</b>	Περιφέρεια βραχίονα στη μασχάλη	Περιφέρεια του αντιβράχιου (στο παχύτερο σημείο)	Περιφέρεια καρπού
<b>cs</b>	Obvod v podpaží	Obvod v předlokti (v nejširším místě)	Obvod zápěstí
<b>lv</b>	Augšdelma augšdaļas apkārtmērs	Apakšdelma apkārtmērs (platākajā punktā)	Plaukstas locītavas apkārtmērs
<b>lt</b>	Rankos apimtis ties pažastimi	Žasto apimtis (storiausioje vietoje)	Riešo apimtis
<b>et</b>	Õlavarre übermõõt kaenla alt	Küünarvarre übermõõt (kõige jämedamast kohast)	Randme übermõõt
<b>sl</b>	Obseg pazduhe	Obseg podlakti (na najširšem delu)	Obseg zapestja
<b>sk</b>	Axilárny obvod ramena	Obvod predlaktia (najširší bod)	Obvod zápästia
<b>hu</b>	A kar kerülete a hónaljnyalban mérve	Az alkar kerülete (a legszélesebb részen)	A csukló kerülete
<b>bg</b>	Обиколка на мишницата	Обиколка на предмишницата (в най-широката част)	Обиколка на китка
<b>ro</b>	Circonfereința brațului la axilă	Circonfereința antebrățului (zona cel mai groasă)	Circonfereința încheieturii
<b>hr</b>	Opseg ruke ispod pazuha	Opseg ruke u laktu (najveći)	Opseg ručnog zgloba
<b>zh</b>	腋下臂围	前臂围 (最粗处)	腕围
<b>ja</b>	脇下の腕周り	前腕周り (いちばん太い部分)	手首周り
<b>ko</b>	겨드랑이 팔 둘레	팔뚝 둘레(가장 굵은 부분)	손목 둘레
<b>ar</b>	محيط الذراع الإبطي	محيط الساعد (في أعرض مكان)	محيط المعصم

Hand - Hand - Hand - Dłoń - Main - Hand - Mano - Mano - Mão - Hånd - Käsi - Hånd -  
 Άκρο χέρι - Ruka - Plauksta - Plaštaka - Kaelaba - Dian - Ruka - Kéz - Рѹка (китка до  
 пръсти) - Mână - Šaka - 手部 - 手 - 손 - اليد

	1	14 - 19 cm
		5.5 - 7.5"
	2	19 - 24 cm
		7.5 - 9.5"
	3	24 - 29 cm
		9.5 - 11.5"

	cA
<b>en</b>	Hand circumference
<b>nl</b>	Omtrek van de hand
<b>sv</b>	Omkrets hand
<b>pl</b>	Obwód dłoni
<b>fr</b>	Circonférence de main
<b>de</b>	Handumfang
<b>it</b>	Circonferenza mano
<b>es</b>	Contorno de la mano
<b>pt</b>	Circunferência da mão
<b>da</b>	Håndens omkreds
<b>fi</b>	Käden ympärysmitta
<b>no</b>	Håndens omkrets
<b>el</b>	Περιφέρεια άκρου χεριού
<b>cs</b>	Obvod ruky
<b>lv</b>	Plaukstas apkārtmērs
<b>lt</b>	Plaštakos apimtis
<b>et</b>	Käelaba ümbermõõt
<b>sl</b>	Obseg dlani
<b>sk</b>	Obvod ruky
<b>hu</b>	A kéz kerülete
<b>bg</b>	Обиколка на ръката
<b>ro</b>	Circumferința mâinii
<b>hr</b>	Opseg šake
<b>zh</b>	手圍
<b>ja</b>	手のひら周り
<b>ko</b>	손 둘레
<b>ar</b>	محيط اليد

## en ADJUSTABLE COMPRESSION DEVICE FOR ARM AND HAND

### Description/Destination

The device comprises a number of adjustable and repositionable overlapping straps which are closed with hook-and-loop fasteners.

Different models are available to fit various parts of the body: foot, calf, knee, thigh, arm and hand.

The calf and thigh models have a tightening system used to control and reproduce the desired level of pressure: low (minimum 20 mmHg), medium (minimum 30 mmHg) and high (minimum 40 mmHg).

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

### Composition

Components for the arm model: polyamide - elastane - polyurethane - polyester.

Components for the hand model: polyamide - synthetic rubber (CR) - polyester.

### Properties/Mode of action

The devices provide compression by applying pressure to the limb. They promote venous and lymphatic circulation, help prevent the formation of edema and maintain volume. Inelastic short-stretch material (average pressure between 20 and 50 mmHg).

### Indications

Management of edema (venous, lymphatic and caused by injury).

Symptomatic treatment of lipedema.

### Contraindications

Do not use in case of major dermatological conditions on an affected limb.

Do not use in the event of known allergy to any of the components.

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) <0.6.

Do not use in the event of decompensated heart failure.

Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg).

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb.

Do not use in case of extra-anatomical bypass of the affected limb.

Do not use in the event of septic thrombosis.

Do not use on the upper limbs in the following cases:

- Do not use in case of brachial plexus condition.
- Do not use in case of vasculitis of the extremities.

### Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

If the device seems to be too tight at the recommended tightness, remove it and see a healthcare professional.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional.

Examine the shape of the leg; if there are any deformed or bony areas, use the protective foam bandages or Varico pads to make the shape of the limb uniform. In the event of a venous ulcer, a Doppler examination is necessary before using the device in order to ensure there is no associated severe arterial disease.

In the event of lipedema, compression can be used to help relieve the symptoms. Seek advice from a healthcare professional for compression and other solutions for managing lipedema.

For some conditions (or situations) such as venous thrombosis, the device should be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product.

Do not apply the product directly to an open wound without a dressing.

### Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness.

Any serious device-related incidents occurring should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

### Instructions for use/Application

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

### Putting on the arm model:

1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up.Ⓐ
2. Slip the arm into the sleeve, Thuasne® logo upwards, sleeve on the inside of the arm.Ⓐ
3. Position the bottom strap Ⓐ: hold the inner part against the skin, tighten the other part over it and fasten with the hook-and-loop fastener (do not fasten to the sleeve) Ⓐ
4. Unroll the other straps (Ⓑ, Ⓒ, Ⓓ) from bottom to top.Ⓒ
5. If necessary, wrap the removable strap Ⓔ over the elbow.Ⓒ  
The hook-and-loop fastener must not be closed over the sleeve.
6. To remove, undo the straps from top to bottom, roll the hook-and-loop fasteners up again and pull the product downwards.

### When used with Mobiderm®:

- For the arm, it is advisable to use Mobiderm® big blocks pads rather than bandages.
- Protect the skin with a suitable protective liner (bandage or other).
- Wrap the Mobiderm® pad around the limb.
- Position the Biflex Self Adjust device on top.
- Adding a Mobiderm® big block pad increases the circumference by an average of 8 cm. This increase must be taken into account when choosing the size of the device.

### Putting the hand model on:

1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up.Ⓐ
2. Open the product.Ⓐ
3. Slide the thumb into the thumb hole (right hand/left hand) Ⓐ
4. Close the hook-and-loop fastener Ⓐ Ⓒ
5. Unroll the Ⓒ bottom strap.Ⓒ
6. Unroll the Ⓒ top strap.Ⓒ
7. To remove, undo the straps from top to bottom and roll the hook-and-loop fasteners up again.Ⓒ

When putting it on and wearing it, ensure no large folds appear.

If necessary, the segment pressure can be adjusted by adjusting the loops from bottom to top.

### When used with Mobiderm®:

- For the extremities (hands/feet), Mobiderm® small block pads should be used.
- Protect the skin with a suitable protective liner (bandage or other).
- Cut the Mobiderm® small blocks pad to the desired size and position it over the protective liner at the site of the edema.
- Position the Biflex Self Adjust device on top.

When putting it on and wearing it, ensure no large folds appear.

If necessary, the segment pressure can be adjusted by adjusting the loops from bottom to top.

### Care

Compression devices are designed for daily use over a 6-month period (i.e. 90 washes). After this time, garments should be replaced to ensure the treatment is effective and technical performance is not affected. Close the hook-and-loop fasteners before washing. Machine washable at 30°C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

### Storage

Store at room temperature, preferably in the original packaging.

### Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2058701 (2025-07)
**THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186