






THUASNE®

Biflex Self Adjust

en	Adjustable compression device for thigh and knee	6
nl	Verstelbaar compressiehulpmiddel voor dij en knie	6
sv	Justerbar kompressionsanordning för lår och knä	7
pl	Regulowana opaska uciskowa na udo lub kolano	8
fr	Dispositif compressif ajustable pour cuisse et genou	9
de	Verstellbare Kompressionsbandage für Oberschenkel und Knie	10
it	Dispositivo di compressione regolabile per coscia e ginocchio	11
es	Dispositivo compresivo ajustable para muslo y rodilla	12
pt	Dispositivo de compressão ajustável para a coxa e o joelho	13
da	Justerbar kompressionsanordning til lår og knæ	14
fi	Säädettävä kompressiotuote reiteen ja polveen	15
no	Justerbar kompresjonsinnretning for lår og kne	16
el	Ρυθμιζόμενο σύστημα συμπίεσης μηρού και γονάτου	17
cs	Nastavitelná kompresní pomůcka pro stehno a koleno	18
lv	Pielāgojams kompresijas izstrādājums augšstilbam un ceļgalam	19
lt	Reguliuojama kompresinė priemonė šlauniai ir keliui	20
et	Reguleeritav kompressiooniseade põlvele	21
sl	Nastavljivi kompresijski pripomoček za stegno in koleno	22
sk	Nastavitelná kompresívna pomôcka na stehno a koleno	22
hu	Állítható kompressziós eszköz a comba és a térdre	23
bg	Регулируемо компресиращо изделие за бедро и коляно	24
ro	Dispozitiv de compresie ajustabil pentru coapsă sau genunchi	25
hr	Podesiva oprema za kompresiju bedra i koljena	26
zh	可调节的大腿或膝盖压力装置	27
ja	大腿部と膝用の調節可能な圧迫装置	28
ko	허벅지와 무릎을 위한 조절 가능한 압박 장치	29
ar	تجهيرة ضاغطة قابلة للتعديل للفخذ والركبة	31







Thigh - Dij - Lår - Udo - Cuisse - Oberschenkel - Coscia - Muslo - Coxa - Lår - Reisi - Lår -
 Μηρός - Stehno - Augšstilbs - Šlaunis - Reis - Stegno - Stehno - Comb - Бедро - Coapsă -
 Bedro - 大腿 - 大腿部 - 허벅지 - ذخفلا


		cE	cG	N	L
	1	34 - 45 cm	46 - 57 cm	$l < 36 \text{ cm} / 14''$	$l > 36 \text{ cm} / 14''$
		13.5 - 18"	18 - 22.5"		
	2	41 - 53 cm	53 - 65 cm		
		16 - 21"	21 - 25.5"		
	3	49 - 60 cm	61 - 72 cm		
		19.5 - 23.5"	24 - 28.5"		
	4	56 - 68 cm	68 - 80 cm		
		22 - 27"	26.5 - 31.5"		

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en ADJUSTABLE COMPRESSION DEVICE FOR THIGH AND KNEE

Description/Destination

The device comprises a number of adjustable and repositionable overlapping straps which are closed with hook-and-loop fasteners.

Different models are available to fit various parts of the body: foot, calf, knee, thigh, arm and hand.

The calf and thigh models have a tightening system used to control and reproduce the desired level of pressure: low (minimum 20 mmHg), medium (minimum 30 mmHg) and high (minimum 40 mmHg).

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Components for the thigh model:

Textile components: polyamide - elastane - polyurethane - polyester.

Tightening system: polyamide - polyurethane.

Components for the knee model:

Polyamide - synthetic rubber (CR) - polyester.

Properties/Mode of action

The devices provide compression by applying pressure to the limb. They promote venous and lymphatic circulation, help prevent the formation of edema and maintain volume.

Inelastic short-stretch material (average pressure between 20 and 50 mmHg).

Indications

Management of edema (venous, lymphatic and caused by injury).

Chronic venous insufficiency (C3 to C6).

Symptomatic treatment of lipedema.

Contraindications

Do not use in case of major dermatological conditions on an affected limb.

Do not use in the event of known allergy to any of the components.

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6.

Do not use in the event of decompensated heart failure.

Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg).

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb.

Do not use in case of extra-anatomical bypass of the affected limb.

Do not use in the event of septic thrombosis.

Do not use on the upper limbs in the following cases:

- Do not use in case of brachial plexus condition.
- Do not use in case of vasculitis of the extremities.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

If the device seems to be too tight at the recommended tightness, remove it and see a healthcare professional.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional.

Examine the shape of the leg: if there are any deformed or bony areas, use the protective foam bandages or Varico pads to make the shape of the limb uniform.

In the event of a venous ulcer, a Doppler examination is necessary before using the device in order to ensure there is no associated severe arterial disease.

In the event of lipedema, compression can be used to help relieve the symptoms. Seek advice from a healthcare professional for compression and other solutions for managing lipedema.

For some conditions (or situations) such as venous thrombosis, the device should be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product.

Do not apply the product directly to an open wound without a dressing.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness.

Any serious device-related incidents occurring should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Putting the thigh model on:

1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up. ①
2. Pull the sleeve ⑥ up to mid-thigh height. ②
3. Unroll the bottom strap ④, tighten the hook-and-loop fastener on the other part of the strap (do not fasten to the sleeve) and tighten on the limb according to the desired level of pressure. ③④
- Take care not to pull on the adjuster.
4. Repeat with all the other straps from bottom to top (⑥, ③ and ④). ④
5. To remove, undo the straps from top to bottom and roll the hook-and-loop fasteners up again. ③

When used with Mobiderm®:

- For the lower limbs, Mobiderm® big blocks pads or bandages should be used.
- Protect the skin with a suitable protective liner (bandage or other).
- Wrap the Mobiderm® pad around the limb or apply the Mobiderm® bandage in a circular fashion around the limb.
- Position the Biflex Self Adjust device on top.
- Adding a Mobiderm® big blocks pad or bandage increases the circumference by an average of 8 cm. This increase must be taken into account when choosing the size of the device.

Using the pressure adjuster ①:

The adjuster comes with 3 levels of pressure: **L** low (light blue line), **M** medium (blue line) and **H** high (dark blue line).

Pressure is reached by pulling on the strap which slides the adjuster until it the outside of the loop is over the desired level of pressure.

When putting it on and wearing it, ensure no large folds appear.

If necessary, the segment pressure can be adjusted by adjusting the loops from bottom to top.

Putting the knee model on:

1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up. ①
2. Position the product on the knee. ③
3. Unroll the middle strap ⑤ and fasten on the other side while tightening it on the limb ④
4. Position the bottom strap ⑥: hold the inner part against the skin, tighten the other part over it and fasten with the hook-and-loop fastener. ③
5. Repeat with the top strap ②. ①
6. To remove, undo the straps from top to bottom and roll the hook-and-loop fasteners up again. ③

When used with Mobiderm®:

- For the lower limbs, Mobiderm® big blocks pads or bandages should be used.
- Protect the skin with a suitable protective liner (bandage or other).
- Wrap the Mobiderm® pad around the limb or apply the Mobiderm® bandage in a circular fashion around the limb.
- Position the Biflex Self Adjust device on top.
- Adding a Mobiderm® big blocks pad or bandage increases the circumference by an average of 8 cm. This increase must be taken into account when choosing the size of the device.

Care

Compression devices are designed for daily use over a 6-month period (i.e. 90 washes). After this time, garments should be replaced to ensure the treatment is effective and technical performance is not affected. Close the hook-and-loop fasteners before washing. Machine washable at 30°C (delicate program). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2058801 (2025-07)
THUASNE UK Ltd
Thuasne SAS - SIREN/RCS Nanterre 542 091 186