



**THUASNE®**

## SpryStep® One

<b>fr</b> Releveur de pied postérieur dynamique.....	3
<b>en</b> Posterior dynamic ankle-foot orthosis (AFO).....	3
<b>de</b> Dynamische Peronäusorthese.....	4
<b>nl</b> Dynamische enkel-voetorthese posterior.....	5
<b>it</b> Tutore dinamico posteriore.....	6
<b>es</b> Antiequino posterior dinámico.....	6
<b>pt</b> Suporte de pé posterior dinâmico.....	7
<b>da</b> Dynamisk posterior ankel-fod-skinne.....	8
<b>fi</b> Dynaaminen posteriorinen jalannostin.....	8
<b>sv</b> Dynamisk posterior droppfotsortos.....	9
<b>el</b> Δυναμικός οπίσθιος ανυψωτήρας ποδιού.....	10
<b>cs</b> Dynamická posteriorní ortéza.....	11
<b>pl</b> Dynamiczna orteza tylna do podtrzymywania stopy.....	11
<b>lv</b> Aizmugurējā stiprinājuma dinamiskā krītošās pēdas ortoze.....	12
<b>lt</b> Dinaminis kulninis pėdos įtvaras.....	13
<b>et</b> Põia dünaamiline posterioorne tugi.....	13
<b>sl</b> Dinamična opornica za dvig stopala.....	14
<b>sk</b> Dynamická zadná ortéza na zdvíhanie chodidla.....	15
<b>hu</b> Hátralról támasztó dinamikus boka-lábfej ortézis.....	15
<b>bg</b> Динамична ортеза за задна стабилизация на стъпалото.....	16
<b>ro</b> Ridicător dinamic posterior pentru picior.....	17
<b>ru</b> Динамический ортез-стоподержатель с шиной на задней поверхности голени.....	17
<b>hr</b> Ortoza za dinamično podizanje stražnjeg dijela stopala.....	18
<b>zh</b> 后护壳式动态踝足抬脚器.....	19
<b>ar</b> رافع نشط لمؤخرة القدم.....	20

<b>fr</b>	Stabilisation	Correction biomécanique	Restitution d'énergie
<b>en</b>	Stabilization	Biomechanical correction	Energy return
<b>de</b>	Stabilisierung	Biomechanische Korrektur	Energierückgabe
<b>nl</b>	Stabilisatie	Biomechanische correctie	Energie-teruggave
<b>it</b>	Stabilizzazione	Correzione biomeccanica	Restituzione di energia
<b>es</b>	Estabilización	Corrección biomecánica	Restitución de energía
<b>pt</b>	Estabilização	Correção biomecânica	Restituição de energia
<b>da</b>	Stabilisering	Biomekanisk korrektion	Genoprettelse af energi
<b>fi</b>	Stabiloi	Biomekaaninen korjaus	Energianpalautus
<b>sv</b>	Stabilisering	Biomekanisk korrigering	Frigör energi
<b>el</b>	Σταθεροποίηση	Βιο-μηχανικός έλεγχος	Αποκατάσταση ενέργειας
<b>cs</b>	Stabilizace	Biomechanická korekce	Restituce energie
<b>pl</b>	Stabilizacja	Korekcja biomechaniczna	Odzyskanie energii
<b>lv</b>	Stabilizēšana	Biomehāniskā korekcija	Energijas atgūšana
<b>lt</b>	Stabilizavimas	Mechaninė biologinė korekcija	Energijos atgavimas
<b>et</b>	Stabiliseerimine	Biomehaaniline korrektsioon	Energia taastamine
<b>sl</b>	Stabiliziranje	Popravek biomehanske motnje	Obnovitev energije
<b>sk</b>	Stabilizácia	Biomechanická korekcia	Obnova energie
<b>hu</b>	Stabilizálás	Biomechanikus korrekció	Energia visszanyerése
<b>bg</b>	Стабилизация	Биомеханична корекция	Възвръщане на енергията
<b>ro</b>	Stabilizare	Corecție biomecanică	Recuperarea energiei
<b>ru</b>	Стабилизация	Биомеханическая коррекция	Восстановление энергии
<b>hr</b>	Stabilizacija	Biomehanička korekcija	Obnova energije
<b>zh</b>	稳定	生物力学矫正	体能恢复
<b>ar</b>	الاستقرار	تصحيح ميكانيكي حيوي	استرجاع الطاقة

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### Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table. This device is a posterior dynamic ankle-foot orthosis (AFO) that supports and/or stabilizes the foot and ankle while walking.

### Composition

Rigid components: carbon fibre - epoxy resin - high density polyethylene. Textile components: polyamide - elastane - polyurethane - ethylene vinyl acetate.

### Properties/Mode of action

The AFO is composed of two parts (rigid **(a)** and soft **(b)** part) already assembled. The rigid part **(a)**, made of composite materials, is positioned under the foot and along the leg to provide stabilisation, biomechanical correction and energy return.

The blue translucent area on the footplate and the shell is trimmable **(b)**. The textile components (foam pads and strap(s)) protect the leg and allow correct position of the device on the limb **(c)**.

## Indications

These indications are biomechanical deficits of neurological, traumatic or muscular origin.  
Weak dorsiflexors.  
Fatiguable footdrop.  
Footdrop.  
Footslap.  
Plantar flexor strength s4.  
Mild spasticity of the foot and ankle.  
Mild knee instability during stance phase.  
Knee hyperextension.  
Circumduction.  
Vaulting (plantarflexion of the contralateral ankle joint).  
High knee gait.  
Hip hiking/contralateral trunk bending.

## Contraindications

Do not use the product if the diagnosis has not been confirmed.  
Do not apply the product in direct contact with broken skin.  
Do not use in the event of known allergy to any of the components.  
Do not use for patients weighing > 135 kg / 300 lbs.  
Severe loss of sensation in the lower limb.  
Open ulcers of the foot, ankle or lower leg.  
Moderate to severe edema of the affected limb.  
Moderate to severe spasticity of the foot and ankle.  
Triplanar instability.  
Plantarflexion contracture.  
Moderate to severe foot deformities.  
Moderate to severe ankle instabilities.  
Running/high impact activities.

## Precautions

Verify the product's integrity before every use.  
Do not use the device if it is damaged.  
The initial fitting and adjustment must be done by a healthcare professional.  
Strictly comply with your healthcare professional's prescription and recommendations for use.  
Check the condition of the affected limb and the state of the skin daily (with particular attention for patients with sensory deficit).  
In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.  
For hygiene, security and performance reasons, do not re-use the product for another patient.  
Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...)  
Do not wear the device in a medical imaging machine.  
The ability to drive a vehicle with the device must be assessed by a healthcare professional and according to local regulations.  
The systematic use of a sock is recommended when wearing the device.  
It is recommended to adequately tighten the device to achieve a good fit on the limb without restricting blood circulation.  
Do not expose the product to extreme temperatures.  
Do not kneel or squat with the device. Ⓞ  
Avoid excessive pressure on the forefoot area:  
- Always place the entire foot on any step or uneven surface. Ⓞ  
- Transition sitting/standing position (chair, toilet, car, ...): put the foot flat on the ground before moving to standing. Use any fixed support (armrests, support bar...) to limit the overload of the foot lifter. Ⓞ

## Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.  
Possible risk of venous thrombosis.  
Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

## Instructions for use/Application

The device must be worn in shoes with the following features:  
- Stiff posterior heel cup sufficiently high to fully encapsulate the foot and the orthosis: recommended height 7 - 7,5 cm / 2 ¼ - 3" Ⓞ;  
- Low heel height (< 2 cm / ¾", 1,1 - 1,3 cm / ¼ - ½" higher than the forefoot area) Ⓞ;  
- Shoe closure system: self-fastening tabs or laces.  
Trainers or athletic shoes are the best type of shoes to use with the device.  
Preparation of the AFO performed by the healthcare professional:  
Ensure that the model (right or left) is appropriate for the patient's needs. Make sure that the patient's footwear complies with the recommendations. The healthcare professional must supervise the fitting of the product and the specific walking conditions of the patient when using the device for the first time.  
Sole pattern and trimming (if necessary):  
Note: only the blue soft area at the tip of the foot, from line to the end of the footplate, can be trimmed.  
Measure the foot length from the heel to the longest toe and/or ask the patient's shoe size.  
If the shoe has a removable inner sole, use it to draw the cutting outline on the device's sole.  
If there is no removable inner sole, position the patient's foot to draw the cutting outline on the device's sole. ①  
Use scissors to cut out the device sole along the traced outline. ②  
Use a buffing cone or fine grade sandpaper to create a smooth finish on trimmed area(s).

Shell pattern and trimming (if necessary):  
Note: only the blue soft areas at the top and at the bottom of the shell, can be cut at the level of the lines.  
Measure the height from the floor to the fibular head.  
Use scissors to cut out the shell at the level of one of the two lines to ensure that the top of the device is at least 2 cm / ¾" below the fibular head. ③  
The lower part of the shell can be cut independently medially and laterally using scissors at the level of one of the two lines.  
Ensure that the medial and lateral heights of the shell correspond to the patient's morphology (height and possible varus or valgus to be corrected).  
Use a buffing cone or fine grade sandpaper to create a smooth finish on trimmed area(s).

Self-fastening tabs placement:  
Stick the 3 rectangular self-fastening tabs to the inside of the shell: the longest in the middle and the other two at the ends. ④  
Stick the 2 round self-fastening tabs to the outside of the shell at these two ends. ⑤  
Shell padding trimming (if necessary):  
Cut the shell padding at the grooves according to the cutting configuration made on the shell. ⑥  
Position the shell padding inside the device shell.  
Fitting the AFO:  
Take the inlay out of the shoe, if removable.  
Place the device into the shoe. ⑦  
Ensure the device's foot plate sits flat on the sole of the shoe and that the posterior heel cup of the shoe is not distorted. ⑧  
Replace the removable inlay on top of the device's footplate, unless it makes the shoe too tight. ⑨  
If there is no removable inlay then disregard this step.  
Loosen the shoe laces and slide the foot inside the shoe. ⑩

You may use a shoe horn if required.

## Check the fit:

Secure the strap(s):  
- Position the buckle on the medial side of the leg.  
- Feed the strap(s) through the corresponding buckle(s). Ⓞ  
- Position the self-fastening tab on the extremity of the strap without the buckle.  
- Ensure that the strap padding is correctly positioned on the front of the tibia.  
- If a strap is too long, remove the self-fastening tab, trim the strap with scissors and reposition the self-fastening tab.  
- Tighten the strap(s).

Ensure laces or other tightening mechanisms are firmly fastened.  
Ensure comfort of foot and leg, with no impingements, prior to use. ⑪  
Additional accessories/spare parts are available to order.

## Fitting the spare parts (by a healthcare professional or by the patient):

The spare parts kit contains the following components : foam pad(s), strap(s), self-fastening tab(s).  
Cut the new shell padding at the grooves according to the cutting configuration made on the shell. ⑫  
Remove the textile parts and the self-fastening tab(s) (if damaged) stuck on the rigid part.  
Clean the area where the self-fastening tabs were applied.  
Replace the self-fastening tab by new one(s) and then position the new foam pad.  
If needed, shorten the strap(s): remove the self-fastening tab, trim the strap(s) and replace the tab.  
Position the buckle of the strap on the rigid part and follow the fitting instructions.

## Care

Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label.  
If the device comes into contact with water, dry the textile part and wipe the rigid part well with a dry cloth.  
If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.  
Rigid components:  
Wash the rigid part with a moist cloth.  
Textile components:  
The soft part can be fully removed for washing.  
Replace in the original location before next use.  
Machine washable at 30°C (delicate programme). Remove the self-fastening tabs before washing. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not dry clean. Do not tumble-dry. Do not iron. Squeeze out excess water. Dry flat. Dry away from any direct heat source (radiator, sun, etc.).

## Storage

Store at room temperature, preferably in the original packaging.

## Disposal

Dispose of in accordance with local regulations.

## COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS

Thuasne® offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:  
- six months for the textile components;  
- two years for the rigid components.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:  
- any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,  
- damages occurred as part of attempts to modify the product.  
Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne® entity.

Any warranty claim will first be reviewed by Thuasne® to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.  
To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.  
If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

## FOR US ONLY

The following guide demonstrates steps to take when minimal self-adjustment is required at time of delivery (as defined by Code of Federal Regulations) HCPCs code L1952 should be used. When more than minimal self-adjustment by a person of expertise at the time of delivery is executed, with appropriate supporting documentation, then HCPCs code L1951 should be used. See current LCD policy for most up to date definitions.

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©Thuasne® - 2048203 (2025-10)
THUASNE UK Ltd