

THUASNE®

MOBIDERM®

fr	Plaque et bandage mobilisateurs pour les œdèmes	3
en	Mobilizing pad and bandage for edema	4
de	Mobilisierende Pads und Bandagen zur Behandlung von Ödemen	6
nl	Mobilisatiepad en -verband voor oedemen	8
it	Placca e benda di mobilizzazione per edemi	9
es	Placa y venda movilizadoras para edemas	11
pt	Placa e ligadura de mobilização para os edemas	13
da	Mobiliserende plade og bind til ødemer	14
fi	Mobilisoiva levy ja sidos turvotukseen	16
no	Kompresjonsbandasjer for ødemer	18
sv	Mobiliserande ark och bandage för ödem	19
el	Πλάκα και επίδεσμος κινητοποίησης για τα οίδηματα	21
cs	Mobilizační podložka a bandáž pro léčbu edémů	23
pl	Arkusz i bandaż mobilizujący na obrzęki	24
lt	Poodinė cirkuliaciją skatinanti plokštelė ir juosta ištinusioms vietoms	26
et	Nahaaluseid kudesid mobiliseeriv padi ja side ödeemide puhuks	27
sl	Mobilizacijska blazinica in povoj za edeme	29
sk	Podložka a obväz na mobilizáciu edémov	31
hu	Ödémacsökkentő mobilizáló lap és pólya	32
bg	Мобилизиращи подложки и бинтове за отоци	34
ro	Placă și bandaj de mobilizare pentru edeme	36
ru	Мобилизующие компрессы и повязки для лечения отеков	37
hr	Podložak i zavoj za imobilizaciju kod edema	39
zh	用于水肿的活动泡沫片和活动绷带	41
ja	リンパ浮腫用活性化プレートおよび包帯	42
ko	부종 모빌라이징 패드 및 붕대	43
ar	رقادة وضادة لتخفيف الوذمات	47

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Bruksanvisning - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - 使用方法 - 사용설명서 - إرشادات الاستعمال

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Conserver cette notice.

en

MOBILIZING PAD AND BANDAGE FOR EDEMA

Description/Use

Mobiderm® mobilizing devices are intended for the treatment of edema of the upper or lower limbs.

Target population: adult and children patients suffering from the indicated pathologies.

This range is available in several versions:

Mobiderm® bandages:

- small blocks (5 x 5 mm) **A**
- big blocks (15 x 15 mm) **B**

Mobiderm® pads:

- small blocks (5 x 5 mm) **A**
- big blocks (15 x 15 mm) **B**

The devices feature Mobiderm® technology, which consists of foam blocks encased between two sheets of non-woven material.

Mobiderm® devices are available in different sizes of pads and bandages and in two foam block sizes to be used according to the needs and the body area treated.

Composition

Textile components: polyester - acrylic adhesive mass.

Foam components: polyurethane (large blocks) - polyethylene (small blocks).

Properties/Mode of action

Mobiderm® devices are intended to be used under a compression bandage.

They provide subcutaneous drainage via Mobiderm® technology, which consists of foam blocks randomly positioned between two sheets of non-woven material.

In contact with the skin, the blocks provide a pressure differential between the contact area beneath them and their periphery, creating multi-directional circulation channels for the lymph under the skin.

The lymph travels from the areas pressed by the blocks towards the free areas.

Efficient lymph drainage and circulation are thus promoted. **C**

Indications

Acute or chronic edema management (e.g., lymphedema).

Edema linked to venous ulcer (CEAP classification stage C6).

Contraindications

Do not use in the event of known allergy to any of the components.

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of septic thrombosis.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of the affected limb.

Specific to the lower limbs:

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6.

Do not use in the event of decompensated heart failure.

Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg).

Do not use in case of extra-anatomical bypass of the affected limb.

Specific to the upper limbs:

Do not use in case of a brachial plexus pathology.

Do not use in case of vasculitis of the extremities.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

It is recommended that a healthcare professional supervises the first application.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

In case of major dermatological conditions affecting a limb, consult a healthcare professional before use.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product,

remove the product and seek the advice of a healthcare professional.

Do not apply the product directly to an open wound without a dressing.

Application of this device can cause the appearance of marks on the skin. These marks are benign and temporary; do not stop your treatment without medical advice.

Due to the constant evolution of the pathology, regular follow-up with a healthcare professional is recommended to check that the compression device is still suitable for the pathology.

Do not use when swimming.

Do not use products on the skin (creams, ointments, oils, gels, patches, etc.) before applying the device as this could damage the product.

In the event of associated advanced venous disease or venous ulcer, seek the advice of a healthcare professional for appropriate compression treatment on the entire leg.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness.

Any serious device-related incidents occurring should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that your device be washed before it is used for the first time.

Verify the product's integrity before every use.

A cotton bandage or protective bandage must be placed underneath the Mobiderm® device.

It is possible to place cushioning devices (such as Varico, foam bandage) underneath Mobiderm® devices.

It is possible to cut the devices to fit the area to be treated according to needs.

It is recommended that the devices be cut with scissors between 2 rows of foam blocks.

Check that there are no folds in the device once in place.

Application of Mobiderm® bandages:

Apply the bandage circularly along the length of the limb without pressure and with an overlap of one row of blocks between each turn.

On the joints, we recommend an overlap of 2 rows of foam blocks in order to maintain the position of the device during various body movements.

Cut off surplus bandage if necessary.

Then apply the compression bandage according to the manufacturer's recommendations.

Application of Mobiderm® pads:

Locate the area to be treated then cut the pad according to the shape of the chosen surface. Apply the pad to the area to be treated without pressure.

Then apply the compression bandage according to the manufacturer's recommendations.

Care

To ensure that Mobiderm® devices last as long as possible and to prevent infections, it is essential to wash them every day, or after each use. Soak the device in tepid soapy water, rub with a clean sponge and squeeze out excess water. Do not use chlorine-based detergents or harsh products. Air dry, away from direct heat (e.g. radiator, sun, etc.). Do not tumble-dry. Do not iron.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2049802 (2025-11)
THUASNE UK Ltd**
Thuasne SAS – SIREN/RCS Nanterre 542 091 186