
THUASNE®

Cicatrex

de Gebrauchsanweisung - Cicatrex Hüfthose
en Instructions for use - Cicatrex hip trouser
fr Consignes d'utilisation - Cuissard compressif Cicatrex
it Istruzioni per l'uso - Pantaloni a vita bassa Cicatrex
es Instrucciones de utilización - Pantalón de compresión de cadera Cicatrex
pt Instruções de uso - Calção de compressão Cicatrex
nl Gebruiksaanwijzing - Cicatrex heupbroek
da Brugsanvisning - Cicatrex hoftebukser
sv Bruksanvisning - Cicatrex höftbyxor
fi Käyttöohje - Cicatrex-lonkkahousut
no Bruksanvisning - Cicatrex hoftebukser
el Οδηγίες χρήσης - Παντελόνι ισχίου Cicatrex
cs Návod k použití - kalhoty Cicatrex
pl Instrukcja użytkowania - spodnie na biodra Cicatrex
et Kasutusjuhend - Cicatrex puusapüksid
lv Lietošanas instrukcija - Cicatrex kompresijas bikses gurnu zonai
lt Naudojimo instrukcija - „Cicatrex“ kompresinės klubų kelnės
sk Návod na použitie - Cicatrex nohavice na bedrá
sl Navodila za uporabo - medicinske hlače za kolke Cicatrex
hu Használati utasítás - Cicatrex csípőszorító
ro Instrucțiuni de utilizare - Pantaloni scurți Cicatrex
bg Инструкция за употреба - Следоперативен компресионен клин Cicatrex (таз и бедра)
hr Upute za uporabu - Cicatrex stezne gaćice
ru Инструкция по использованию - бедренный чулок Cicatrex
tr Kullanım kılavuzu - Cicatrex kalça bandajı
ko 사용 설명서 - Cicatrex 힙 압박 팬츠
ja 取扱説明書 - Cicatrex ヒップパンツ
ar تعليمات الاستخدام - ملابس ضاغطة للثوب

Circumferences in cm	Sizes	S	M	L	XL	2XL	3XL	4XL
Hip		78-90	90-109	109-128	128-135	135-142	142-149	149-156
Thigh (measured about 2 cm below the crotch)		30-42	42-54	54-64	64-69	69-73	73-76	76-81

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34% Viskose, 31% Polyamid, 23% Polyester, 8% Elasthan, 4% Polyurethan.

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Product description

The Cicatrex hip trouser is used post-operatively to provide compression to the operated area. It consists of an adjustable belt, a leg section and a pad to apply compression to the scar area.

Purpose

Compression textiles are basic elements in scar tissue therapy. The objective is to obtain flat and elastic scar tissue which does not reduce the quality of life. Cicatrex compression garment is used for the treatment of scars. The aim is to promote healthy scar maturation and development through the even, planar pressure exerted by the compression garments. It compensates for the natural resilience and surface elasticity of the skin. It prevents the scar from rising above the level of the skin. It reduces blood circulation in the scar. It supports the parallel alignment of collagen fibres. It protects the fresh scar against external influences and relieves itching. It prevents inflammatory reactions in the tissue. It improves the appearance and function of the scar. Compression therapy generally is a long-term treatment and requires the corresponding products which – if the medical indication is given – the patient him/herself applies upon the doctor's instruction.

Indications

Damage of the skin and the skin's healing function with hypertrophic scarring or keloid formation.

- Following injuries / burns / surgical interventions
- Skin grafts
- For improvement of scar / wound healing until complete healing of the scar
- For the stabilisation of surgical results

Absolute contraindications

Peripheral arterial occlusive disease (PAOD) stage III/IV; decompensated heart failure; septic phlebitis; phlegmasia coerulea dolens; acute erysipelas; acute bacterial, viral or allergic inflammations; untreated or advanced dermatoses /contact dermatitis; untreated open wounds; swelling of extremities of unknown origin; untreated malignant tumours.

Relative contraindications

Pronounced exudative dermatoses; malignant lymph oedema; peripheral arterial occlusive disease (PAOD) stage I/II; peripheral neuropathy associated with diabetes mellitus; sensitivity disorders of extremities; compensated heart failure; skin infections; dermatoses (e.g. blistering dermatoses); intolerance of pressure or of product ingredients; primary chronic polyarthritis. In case of open wounds in the treatment area, the wound area must first be covered with an adequate wound dressing or a bandage before the compression garment is applied.

Thuasne Deutschland GmbH accepts no liability for failure to observe these contraindications.

Side effects

No side effects are known if properly used. In case of any negative changes during the application, please consult your doctor or the medical retailer immediately. In case of complaints in connection with the product, such as damage of the knitted fabric or lack of fit, please contact directly your medical retailer.

The manufacturer shall not be liable for damages / injuries which resulted due to improper handling or misuse.

Precautions

Follow the instructions given by your doctor, by an orthopaedic technician or medical expert personnel. If you experience severe pain or discomfort, please see your doctor, an orthopaedic technician or medical expert personnel. Store at room temperature, preferably in the original packaging. To guarantee its effectiveness and for reasons of hygiene, the product may not be passed on to other patients. Any serious incident that occurs in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient are/is domiciled.

Instructions for use

The hip trouser can be used for either the left or right leg. 2 body measurements are required for fitting. The elastic Cicatrex material and the Velcro fastenings allow for adjustment across a wide range of sizes. If the measurements differ significantly, you can use the waist strap (1) from a different size. To open it, the Velcro fastening on the leg section (2) can be fully undone. This Velcro fastening allows for quick inspection of the area operated on. Position the pad (3) in the centre of the scar area. Then secure the leg section (2) using the Velcro fastener.

Care

Compression garments should be washed daily, if possible. The garment fabric not only collects remainders of ointments and grease but also skin scales, sweat and dust. Regular washing will protect you from infections and your garment from wear. You need a spare set of garments so that you need not interrupt your compression therapy while the other set is washed and dried.

Always turn your compression garment inside out to wash it. Use a mild detergent for washing by hand or for machine wash (gentle cycle at 40° C) without optical brightener and fabric softener. In a laundry bag, compression garments are specially protected. In order to reduce the drying period, put the compression garments into a terry cloth towel, roll it up firmly and press the garment out (do not wring it out). Do not leave the moist compression garment in the wet towel. Hang up your garment loosely and leave it to air-dry. Never put it in the sun or onto a radiator to dry.

Do not expose compression garments to any damaging effects, such as: strong sunlight, intense heat, ointments, oils and greases, stain remover, gasoline and aggressive detergents.

Material composition

34% viscose, 31% polyamide, 23% polyester, 8% elastane, 4% polyurethane.

Storage information and duration of use

Store in a dry place and protect from sunlight. For series products, a maximum shelf life of 36 months applies. This is followed by the duration of use of 6 months for the compression garment. The hourglass pictogram shows the duration of use printed on the carton label. Customised manufacture concerns products for immediate use.

For serial sizes and custom-tailored garments, we guarantee their medical efficacy for a maximum period of wear of six months. Prerequisite for this will be the proper treatment of the compression materials (for example, their care or putting them on or taking them off). If necessary, your doctor should then prescribe new compression garments. In case of any significant change of relevant body sizes (e.g. due to a successful therapy, weight change), new measurements and a new supply may also be justified earlier already.

Disposal

You may dispose of your medical compression garment via the residual or household waste. No special criteria exist for their elimination. Please ensure you comply with the applicable local and national regulations.

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Thuasne UK Ltd