THUASNE®

GenuPro Ligastrap®

ADAPT

fr	Genouillère ligamentaire à sanglage fonctionnel avec semi-ouverture et armatures interchangeables.....	8
en	Ligament knee brace with functional strapping, semi-open design, and interchangeable supports.....	10
de	Kniebandage mit funktionellen Zuggurten, halboffenem Design und austauschbaren seitlichen Verstärkungen	12
nl	Ligamentaire kniebrace met functionele banden met halve opening en vervangbare steunen	15
it	Ginocchio per legamenti con tiranti funzionali con semi-apertura e armature intercambiabili. 18	
es	Rodillera ligamentaria con correas funcionales, semiabertura y montantes intercambiables	20
pt	Joelheira ligamentar com correia funcional com semiabertura e reforços intercambiáveis.....	22
da	Ligamentknæbind med funktionel strop, halv åbning og udskiftelige støttestrukturer.....	25
fi	Puoliavoin nivelsidetuki, jossa on toiminnallisuutta edistävät tukisiteet ja keskenään vaihdettavat vahvikkeet	27
sv	Knäskydd för ligament med funktionellat remsystem med halvöppningsbara och utbytbara ramar	30
el	Επιγονατίδα συνδέσμων και επιγονατίδας του γονάτου με λειτουργικό δέσιμο, μισή οπή και εναλλάξιμες ενισχύσεις	32
cs	Kolenní ortéza pro stabilizaci vazů s funkčním zapínáním a polootevřeným designem s vyměnitelnými dlahami	34
pl	Orteza stawu kolanowego na wiązadła z funkcjonalnymi paskami, półotwarcie i wymiennymi ramkami.....	37
lv	Ceļa locītavas saišu ortoze ar funkcionālo saiti, daļēju atveri un nomaināmiem stiprinājumiem.....	39
lt	Kelio raiščių įtvares su funkciniais diržais ir keičiamomis tvirtinimo detalėmis, pusiau atveriamas.....	42
et	Põlveliigese sidemeid toetav funktsionaalse rihmasüsteemiga põlvvetugi, millel on poolavanev kinnis ja vahetatavad tugevdused.....	44
sl	Opornica za kolenske vezi s funkcijskimi trakovi, polodprtino in zamenljivim ogrodjem.....	46
sk	Kolenná ortéza na stabilizáciu väzov s funkčným zapínaním, pootvorením hornej prednej časti a vymeniteľnými dlahami	48
hu	Térdortézis térdszalagok rögzítésére funkcionális pántokkal, félig nyitott kialakítással és cserélhető merevítőekkel	51
bg	Наколенка с функционална лента, с отваряне наполовина и със сменяеми подложки.....	53
ro	Genunchieră ligamentară cu curea funcțională, cu semi-deschidere și cadre interschimbabile.....	56
ru	Наколенник лигаментарный с функциональными ремнями с полуоткрытием и съёмными вставками	58
hr	Ortoza za ligamente koljena s funkcionalnim trakama s poluotvorom i zamjenjivim učvršćivačima	61
zh	带有功能性搭扣的半开放式韧带髌骨护膝，配有互换支具	63
ar	دعامة أربطة الركبة مزودة بحزام وظيفي وفتحة نصفية وإطارات واقية قابلة للتبديل فيما بينها ...	67

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu -

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Conserver cette notice.

en

LIGAMENT KNEE BRACE WITH FUNCTIONAL STRAPPING, SEMI-OPEN DESIGN, AND INTERCHANGEABLE SUPPORTS

Description/Purpose

The device is for the treatment of the indications listed and for patients whose measurements correspond to the sizing table only.

Bilateral product.

Available in 7 sizes.

The device is composed of:

- An anatomically shaped knit fabric with a comfort zone at the back of the popliteal fossa.
- A specific knitted area at the patella (kneecap) ①.
- A functional double strap ②.
- A semi-opening for easy donning and a better fit around the thigh ③.
- A pair of removable rigid articulated supports ④.
- A pair of removable flexible stays ⑤.

This device can be used when resuming sports activities.

Composition

Textile components: polyamide - polyester - polyurethane - elastane - acrylic resin.
Rigid components: steel - thermoplastic elastomer - polypropylene.

Properties/Mode of action

An anatomically-shaped knitting ensuring compression, support and proprioception.

Stabilisation of the knee joint ligaments thanks to double strapping, reinforced by two rigid articulated supports.

Semi-opening on the thigh to make it easier to fit, and so it can be adapted to a wider range of morphologies.

Pull tabs to facilitate application.

Indications

Joint instability/laxity.

Knee pain and/or swelling management.

Only with articulated rigid supports:

Conservative treatment of minor knee ligament injuries (lateral ligaments).

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

If the measurements taken do not correspond to the same size, favor the size corresponding to measurement **(b)**.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient. It is recommended to adequately tighten the device to achieve support without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in color of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

In the event of history of venous or lymphatic disorders, consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

The lateral uprights have been designed for use with the device only. Do not use the knee brace with other uprights.

Do not wear the device in a medical imaging machine.

Do not use the device if certain products are applied to the skin (creams, ointments, oils, gels, patches etc.).

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Fitting

It is recommended to wear the device directly on the skin, unless contraindicated.

The knee brace is supplied with articulated rigid frames already pre-inserted in the slots provided on each side.

Depending on use, leave them in place or replace them with the flexible stays supplied.

To replace the articulated rigid supports with flexible stays:

On the outside of the knee brace, there is an opening in the upper part of the lateral uprights. ❶

Pull down while holding the lower part of the support and the top pull tab to release the upper end of the articulated rigid supports from their covers. ❷

Once removed, insert the flexible stays through the cover openings.

Ensure that the upper end of the covers is properly re-engaged over the flexible stays.

To replace the flexible stays with articulated rigid supports:

Remove the flexible stays and insert the articulated rigid supports into the covers.

Ensure that they are positioned correctly: the shortest part should be positioned facing downwards, bending backwards. ❸

Make sure that the curvature of the ends points outwards so it doesn't press against the leg causing discomfort. ❹

To simplify positioning, identification tabs have been added to the lower external part of the articulated rigid supports.

Fitting the knee brace:

To avoid damaging the knit fabric, ensure that the hook-and-loop fasteners are always attached to the surface designed for this when putting on and taking off the brace.

Loosen the strap so the device is semi-open. A

Slip the knee brace on the leg, with the foot extended, using the pull tabs. B

Ensure correct positioning by checking that the kneecap is centered in relation to the specific knitted area.

Close the strap, adjusting the tightness to the circumference of the thigh. C

Check that the brace is pulled completely taut over the leg.

To assist with fitting, markings are provided on the straps and knee brace. The straps ① have 1 punch hole, the straps ② have 2 punch holes; make sure to position the straps correctly on the corresponding parts of the knee brace.

Pull and position the straps ① such that they pass on either side of the kneecap and attach them to the corresponding hook-and-loop fastener on the knee brace D.

Next, pull and position the straps ② such that they pass on either side of the kneecap and attach them to the corresponding hook-and-loop fastener on the knee brace E.

For comfort, check that there are no creases in the straps.

If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions.

Care

Remove articulated rigid supports before washing, as described in the "Instructions for use" section. Close the self-fastening tabs before washing. Machine washable at 30°C (delicate program). The device has been tested during 20 wash cycles. Use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not dry clean. Do not tumble-dry. Do not iron. Squeeze out any excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, rinse it with clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Para ajudar a colocação, um sistema de marcas está disponível nas correias e na joelheira. As correias © têm 1 marca, as correias ® têm 2 marcas; verificar a posição correta das correias nas partes correspondentes da joelheira.
THUASNE UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186