THUASNE®

SILIMA®

de	**Gebrauchsanweisung - Post-OP BH Cara**
en	**Instructions for use - Post-OP bra Cara**
fr	**Consignes d'utilisation - Soutien-gorge postopératoire Cara**
it	**Istruzioni per l'uso - Reggiseno post-operatorio Cara**
es	**Instrucciones de utilización - Sujetador postoperatorio Cara**
pt	**Instruções de uso - Soutien pós-operatório Cara**
nl	**Gebruiksaanwijzing - Post-OP beha Cara**
da	**Brugsanvisning - Post-OP BH Cara**
sv	**Bruksanvisning - Post-op behå Cara**
fi	**Käyttöohje - Cara - postoperatiiviset kompressiorintaliivit**
no	**Bruksanvisning - Post-OP BH Cara**
el	**Οδηγίες χρήσης - Μετεγχειρητικός στηθόδεσμος Cara**
cs	**Návod k použití - Pooperační podprsenka Cara**
pl	**Sposób użytkowania - Biustonosz pooperacyjny Cara**
et	**Kasutusjuhend - Cara operatsioonijärgne rinnahoidja**
lv	**Lietošanas instrukcija - Pēcoperācijas krūšturis Cara**
lt	**Naudojimo instrukcija - Pooperacinė liemenėlė „Cara“**
sk	**Návod na použitie - Pooperačná podprsenka Cara**
sl	**Navodila za uporabo - Pooperativni nedrček Cara**
hu	**Használati utasítás - Cara operáció utáni melltartó**
ro	**Instrucțiuni de utilizare - Sutien post-operator Cara**
bg	**Инструкция за употреба - Следоперативен сутиен Cara**
hr	**Upute za uporabu - Postoperativni grudnjak Cara**
mk	**Упатства за користење - Post-OP BH Cara**
sr	**Uputstvo za upotrebu - Postoperativni grudnjak Cara**
tr	**Kullanım talimatı: Ameliyat Sonrası Sutyen Cara**
ru	**Инструкция по использованию - Послеоперационный бюстгальтер Cara**
ja	**使用説明書 - 術後ブラジャーCara**
ar	**(Post-OP BH Cara) تعليمات استخدام حمالة الصدر بعد الجراحة**
he	**הוראות שימוש חזיית Cara לאחר ניתוח**
ko	독일어: 사용설명서 **-** Cara 수술 후 브라
zh	使用说明书 **- Post-OP BH Cara**

de
Post-OP BH Cara
Postoperativer Kompressions-BH (Art.-Nr. T57230)

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en

Cara postoperative bra

Postoperative compression bra (Art. No. T57230)

Indication and purpose

The postoperative compression bra is used for uniform compression in the mammary area. Haematomas and postoperative oedemas are minimised. The bra's comfortable handling makes the inspection of wounds easier and reduces pain. Stabilisation and immobilisation of the post-surgical breast is realised after any plastic breast surgery.

The postoperative compression bra consists of a highly elastic and skin-friendly material. The material is doubled in the chest-enclosing front and thus ensures uniform compression.

Triple adjustable straps, hook and eye closures, as well as the front zipper ensure that the bra is easy to adjust and to put on or take off.

The broad, elastic underbust band prevents rolling, pressure and abrasion points, as well as constrictions. The bra is offered in many different sizes. In combination with the postoperative belt (Art. T57231), the bra can be used for the fixation and stabilisation of the breast implant after all plastic breast surgeries.

Indications

Ensuring results after surgical interventions involving:

- Conservative breast surgery
- Breast reconstruction
- Breast reduction
- Breast augmentation

Contraindications

Allergic or inflammatory skin changes (e.g. swelling, reddening of the body area to be cared for). Use only on intact skin, not on infected or acutely inflammatory skin. Please consult with the treating physician regarding the indication and wearing mode of this orthopedic aid.

Sizes

	AA/A B/C D/E	AA/A B/C D/E	AA/A B/C D/E	AA/A B/C D/E	AA/A B/C D/E	AA/A B/C D/E	AA/A B/C D/E	AA/A B/C D/E*	AA/A* B/C*	AA/A* B/C*	AA/A* B/C*
Int/D	70	75	80	85	90	95	100	105	110	115	120
F/E	85	90	95	100	105	110	115	120	125	130	135
I	1	2	3	4	5	6	7	8	9	10	11
GB/US	32	34	36	38	40	42	44	46	48	50	52
AUS/NZ	10	12	14	16	18	20	22	24	26	28	30

*Available only in black colour

How to put on/strap on
The bra is put on and closed in front with a hook and eye closure and the zipper. The length of the straps is triple adjustable.

Precautionary measures
Please follow the instructions by the physician or the orthopedic technician. In case of severe discomfort or unpleasant sensations, please contact your physician or orthopedic technician. Store at room temperature, preferably in its original packaging. Please follow the instructions by the physician or the orthopedic technician. This medical device is intended for the care of an individual female patient. In order to ensure its efficacy and for reasons of hygiene, this device may not be passed on to other female patients. Any serious incidents which occurred in connection with this product are to be reported to the manufacturer and the responsible authority of the member state in which the user and/or the patient are resident.

Washing instructions
60 degrees water temperature; no fabric softener (also see textile label)
Prior to washing, please close all closures!

Warranty
According to the statutory provisions, warranty is extended to material and processing if properly handled according to the instructions for care.

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Thuasne UK Ltd
Thuasne Deutschland GmbH